































1. Select the yoke type; E-200 type yoke for one piece cable with fixed leads, E-201 type yoke for detachable cable with Safety DIN 1.6 mm replaceable leads, E-202 type yoke for detachable cable with dual pole SIEMENS style replaceable leads, E-204 type yoke for detachable cable with dual pole GE style replaceable leads, E-205 type yoke for detachable cable with dual pole HP style replaceable leads, E-206 type yoke for detachable cable with dual pole SPACELABS style replaceable leads and E-207 type yoke for detachable cable with NIHON KOHDEN style replaceable lead group. For PHILIPS NEW type replaceable lead groups select E-208 type yokes.

Table-1 FMT ECG Patient Cable Yokes


 <p>E200-5XXX</p>	
 <p>E200-4XXX</p>	
 <p>E200-3XXX</p>	
 <p>E201-5XXX</p>	
 <p>E201-3XXX</p>	
 <p>E202-5XXX E202N-5XXX*</p>	
 <p>E202-3XXX</p>	
 <p>E204-XXXX</p>	

*IEC-AHA leadwire placement for Datascope cable yokes


Table-1 FMT ECG Patient Cable Yokes (Continued)

<p>▲ E205-5XXX</p>		
<p>▲ E205-3XXX</p>		
<p>▲ E206-5XXX</p>		
<p>▲ E206-3XXX</p>		
<p>▲ E207-3XXX E207-6XXX</p>		
<p>▲ E208-5XXX</p>		
<p>▲ E208-3XXX</p>		















2- Select the patient termination of the leadwires of the cable, “S” for snap, “G” for grabber.













- 3- Select the Color Coding of the leadwires. “A” for AHA color coding, “I” for IEC color coding.
- 4- Select the compatible brand and monitor side connector of the cable by referring to Table-2.
- 5- Complete the part number by putting the desired parameters instead of variables shown by X.




6- Replacement leadwires for E201, 202, 204, 205, 206, 207 and 208 type cables can be bought separately please refer to following chart for part numbers.

▼ Cable Yoke Side	▼ Patient Side	
 ▲ E21/XX		DIN Style Leadwire with Snap Patient End (for E201 yoke cables)
 ▲ E22/XX		DIN Style Leadwire with Grabber Patient End (for E201 yoke cables)
 ▲ E31/XX		SIEMENS/DRAGER Dual pole Leadwire with Snap Patient End (for E202 yoke cables)
 ▲ E32/XX		SIEMENS/DRAGER Dual pole Leadwire with Grabber Patient End (for E202 yoke cables)
 ▲ E41/XX		GE Dual pole Leadwire with Snap Patient End (for E204 yoke cables)
 ▲ E42/XX		GE Dual pole Leadwire with Grabber Patient End (for E204 yoke cables)
 ▲ E51/XX		HP/PHILIPS Dual pole Leadwire with Snap Patient End (for E205 yoke cables)

▼ Cable Yoke Side	▼ Patient Side	
 ▲ E52/XX		HP/PHILIPS Dual pole Leadwire with Grabber Patient End <i>(for E205 yoke cables)</i>
 ▲ E61/XX		SPACELABS Dual pole Leadwire with Snap Patient End <i>(for E206 yoke cables)</i>
 ▲ E62/XX		SPACELABS Dual pole Leadwire with Grabber Patient End <i>(for E206 yoke cables)</i>
 ▲ E91/XX		DATASCOPE Dual pole Leadwire with Snap Patient End
 ▲ E92/XX		DATASCOPE Dual pole Leadwire with Grabber Patient End

Part Numbering of FMT Leadwires: XX is variable for Color Code & Location Designator.
 For IEC Color Code: 0C: White C, 0L: Yellow L, 0R: Red R, 0N: Black N, 0F: Green F. Example: E52/0R
 For AHA Color Code: 0V: Brown V, C LA: Black LA, RA: White RA, RL: Green RL, LL: Red LL. Örnek: E61/RA

▶ Lead Wires for E202, E204 ve E205 Cable Yokes



* E3X/5IN : IEC 5 leadwire setting for Datascope cable yokes.
 ** E3X/3IN : IEC 3 leadwire setting for Datascope cable yokes with 2 dummy leads.































Cable Yoke Side	Patient Side	
	  E71/3I E71/3A	NIHON KOHDEN Leadwire set with Snap Patient End (for E207 yoke cables)
	  E72/3I E72/3A	NIHON KOHDEN Leadwire set with Grabber Patient End (for E207 yoke cables)
	  E71/6I E71/6A	NIHON KOHDEN Leadwire set with Snap Patient End (for E207 yoke cables)
	  E72/6I E72/6A	NIHON KOHDEN Leadwire set with Grabber Patient End (for E207 yoke cables)
	  E81/5I E81/5A	PHILIPS Leadwire set with Snap Patient End (for E208 yoke cables)
	  E82/5I E82/5A	PHILIPS Leadwire set with Grabber Patient End (for E208 yoke cables)
	  E81/3I E81/3A	PHILIPS Leadwire set with Snap Patient End (for E208 yoke cables)
	  E82/3I E82/3A	PHILIPS Leadwire set with Grabber Patient End (for E208 yoke cables)
	  E41/3IT E41/3AT	GE Leadwire set with Snap Patient End (for E204 yoke cables)
	  E42/3IT E42/3AT	GE Leadwire set with Grabber Patient End (for E204 yoke cables)
	  E41/5IT E41/5AT	GE Leadwire set with Snap Patient End (for E204 yoke cables)
	  E42/5IT E42/5AT	GE Leadwire set with Grabber Patient End (for E204 yoke cables)



Table-2 ECG Patient Cables Compatibility Table (Continued)





Connector	Part Number	Yoke Type
 Connector No: 001 BIOCARE		E200-3001/XX One piece (monolithic), 3 Lead
		E201-3001/XX DIN type detachable, 3 Lead
		E202-3001/XX Dual pole detachable, 3 Lead
IM 7E/8E/12E/15E, BM 9000/9000S.		E200-5001/XX One piece (monolithic), 5 Lead
		E201-5001/XX DIN type detachable, 5 Lead
		E202-5001/XX Dual pole detachable, 5 Lead
 Connector No: 051 BIONET		E200-3051/XX One piece (monolithic), 3 Lead
		E201-3051/XX DIN type detachable, 3 Lead
		E202-3051/XX Dual pole detachable, 3 Lead
BM3.		
 Connector No: 140 BIOSYS		E200-3053/XX One piece (monolithic), 3 Lead
		E201-3053/XX DIN type detachable, 3 Lead
		E202-3053/XX Dual pole detachable, 3 Lead
BPM103.		
 Connector No: 022 BRUKER – ODAM		E200-3080/XX One piece (monolithic), 3 Lead
		E201-3080/XX DIN type detachable, 3 Lead
		E202-3080/XX Dual pole detachable, 3 Lead
SM784, SM785, Physiogard 910.		
 Connector No: 058 CORPULS		E200-4058/XX One piece (monolithic), 4 Lead
		E201-4058/XX DIN type detachable, 4 Lead
		E202-4058/XX Dual pole detachable, 4 Lead
Corpuls 08/16.		
 Connector No: 001 CRITICARE SYSTEMS (Csi)		E200-3001/XX One piece (monolithic), 3 Lead
		E201-3001/XX DIN type detachable, 3 Lead
		E202-3001/XX Dual pole detachable, 3 Lead
		E200-5001/XX One piece (monolithic), 5 Lead
		E201-5001/XX DIN type detachable, 5 Lead
		E202-5001/XX Dual pole detachable, 5 Lead
506, 507, 508, 1100, 1508, 504US, 602-4 USP, Poet 2, Scholar.		
 Connector No: 001 DATASCOPE		E200-3001/XX One piece (monolithic), 3 Lead
		E201-3001/XX DIN type detachable, 3 Lead
		E202-3001/XX Dual pole detachable, 3 Lead
		E200-5001/XX One piece (monolithic), 5 Lead
		E201-5001/XX DIN type detachable, 5 Lead
		E202-5001/XX Dual pole detachable, 5 Lead
2000, 2001, 2002, 2100, 2200, 3000,3000A, Baloon Pump, M/D3 w/lead select, M/D3A, M/D4, Passport IL/XG/XL, System 90/90T/97.		
 Connector No: 010 DATEX		E200-3010/XX One piece (monolithic), 3 Lead
		E201-3010/XX DIN type detachable, 3 Lead
		E202-3010/XX Dual pole detachable, 3 Lead
		E200-5010/XX One piece (monolithic), 5 Lead
		E201-5010/XX DIN type detachable, 5 Lead
		E202-5010/XX Dual pole detachable, 5 Lead
AS/3, CS/3, S/5 (with Datex Module), Cardiocap1, Cardiocap 2, Cardiocap 5, Light Monitor, Satelliteplus.		



Table-2 ECG Patient Cables Compatibility Table (Continued)

Connector

Part Number

Yoke Type

 <p>Connector No: 060</p>		E200-3060/XX	One piece (monolithic), 3 Lead
		E201-3060/XX	DIN type detachable, 3 Lead
 <p>Connector No: 001</p>		E202-3060/XX	Dual pole detachable, 3 Lead
		E204-3060/XX	Dual pole GE detachable, 3 Lead
 <p>Connector No: 069</p>		E200-5060/XX	One piece (monolithic), 5 Lead
		E201-5060/XX	DIN type detachable, 5 Lead
 <p>Connector No: 022</p>		E202-5060/XX	Dual pole detachable, 5 Lead
		E204-5060/XX(*)	Dual pole GE detachable, 5 Lead
<p>GE HEALTHCARE Eagle 1000/3000/5000, Solar 8000/9500, Dash 1000/2500/4000/5000, TRAM 100/200/300/400/500 modules, MAC-LAB.</p>		E200-3001/XX	One piece (monolithic), 3 Lead
 <p>Connector No: 001</p>		E201-3001/XX	DIN type detachable, 3 Lead
		E202-3001/XX	Dual pole detachable, 3 Lead
 <p>Connector No: 069</p>		E200-5001/XX	One piece (monolithic), 5 Lead
		E201-5001/XX	DIN type detachable, 5 Lead
 <p>Connector No: 022</p>		E202-5001/XX	Dual pole detachable, 5 Lead
		<p>GOLDWAY UT4000F, UT6000A.</p>	
 <p>Connector No: 022</p>		E200-3069/XX	One piece (monolithic), 3 Lead
		E201-3069/XX	DIN type detachable, 3 Lead
 <p>Connector No: 157</p>		E202-3069/XX	Dual pole detachable, 3 Lead
		E200-5069/XX	One piece (monolithic), 5 Lead
 <p>Connector No: 148</p>		E201-5069/XX	DIN type detachable, 5 Lead
		E202-5069/XX	Dual pole detachable, 5 Lead
<p>HELLIGE Eagle 1000, Cardioserv, SMV and SMV series, Vicom-SM, Servomed4.</p>		E200-3091/XX	One piece (monolithic), 3 Lead
 <p>Connector No: 022</p>		E201-3091/XX	DIN type detachable, 3 Lead
		E202-3091/XX	Dual pole detachable, 3 Lead
<p>HME - HUNTLEIGH Smartsigns Compact 750/1000.</p>		E200-3040/XX	One piece (monolithic), 3 Lead
 <p>Connector No: 157</p>		E201-3040/XX	DIN type detachable, 3 Lead
		E202-3040/XX	Dual pole detachable, 3 Lead
 <p>Connector No: 148</p>		E200-5040/XX	One piece (monolithic), 5 Lead
		E201-5040/XX	DIN type detachable, 5 Lead
 <p>Connector No: 001</p>		E202-5040/XX	Dual pole detachable, 5 Lead
		<p>HME - HUNTLEIGH Innocare-T.</p>	
 <p>Connector No: 001</p>		E200-5006/XX	One piece (monolithic), 5 Lead
		E201-5006/XX	DIN type detachable, 5 Lead
 <p>Connector No: 022</p>		E202-5006/XX	Dual pole detachable, 5 Lead
		<p>INNOMED Millenia, Omni Trak.</p>	
 <p>Connector No: 001</p>		E200-3001/XX	One piece (monolithic), 3 Lead
		E201-3001/XX	DIN type detachable, 3 Lead
 <p>Connector No: 022</p>		E202-3001/XX	Dual pole detachable, 3 Lead
		E200-5001/XX	One piece (monolithic), 5 Lead
 <p>Connector No: 157</p>		E201-5001/XX	DIN type detachable, 5 Lead
		E202-5001/XX	Dual pole detachable, 5 Lead
 <p>Connector No: 022</p>		E200-3022/XX	One piece (monolithic), 3 Lead
		E201-3022/XX	DIN type detachable, 3 Lead
 <p>Connector No: 157</p>		E202-3022/XX	Dual pole detachable, 3 Lead
		E200-5022/XX	One piece (monolithic), 5 Lead
 <p>Connector No: 148</p>		E201-5022/XX	DIN type detachable, 5 Lead
		E202-5022/XX	Dual pole detachable, 5 Lead
<p>KONTRON 7000, 7135, 7135B, 7141, 7210, 7250, 7251, K2000, KAAT, Kolormar, Micromon, Minimon, Superman, Vitalmon.</p>			

(*) E4X and E4X/5X Leadwire sets are available for E204 Cable yokes. Please indicate your selection in brackets after the part number. EXAMPLE: E204-5060/SI (E42/5I)



Table-2 ECG Patient Cables Compatibility Table (Continued)








Connector	Part Number	Yoke Type
 Connector No: 042 KONTRON LM.	E200-3042/XX	One piece (monolithic), 3 Lead
	E201-3042/XX	DIN type detachable, 3 Lead
	E202-3042/XX	Dual pole detachable, 3 Lead
 Connector No: 001 MDE Escort 100/101/102/200/300, Prism.	E200-3001/XX	One piece (monolithic), 3 Lead
	E201-3001/XX	DIN type detachable, 3 Lead
	E202-3001/XX	Dual pole detachable, 3 Lead
	E200-5001/XX	One piece (monolithic), 5 Lead
	E201-5001/XX	DIN type detachable, 5 Lead
E202-5001/XX	Dual pole detachable, 5 Lead	
 Connector No: 243 MEDTRONIC LifePak 20.	E200-3057/XX	One piece (monolithic), 3 Lead
	E200-5057/XX	One piece (monolithic), 5 Lead
 Connector No: 089 MEK MP500/1000.	E200-3089/XX	One piece (monolithic), 3 Lead
	E201-3089/XX	DIN type detachable, 3 Lead
	E202-3089/XX	Dual pole detachable, 3 Lead
	E200-5089/XX	One piece (monolithic), 5 Lead
	E201-5089/XX	DIN type detachable, 5 Lead
E202-5089/XX	Dual pole detachable, 5 Lead	
 Connector No: 080 MEK MP700/800.	E200-3079/XX	One piece (monolithic), 3 Lead
	E201-3079/XX	DIN type detachable, 3 Lead
	E202-3079/XX	Dual pole detachable, 3 Lead
	E200-5079/XX	One piece (monolithic), 5 Lead
	E201-5079/XX	DIN type detachable, 5 Lead
E202-5079/XX	Dual pole detachable, 5 Lead	
 Connector No: 029 MEK MP1000NT, MP1300, MP.	E200-3029/XX	One piece (monolithic), 3 Lead
	E201-3029/XX	DIN type detachable, 3 Lead
	E202-3029/XX	Dual pole detachable, 3 Lead
	E200-5029/XX	One piece (monolithic), 5 Lead
	E201-5029/XX	DIN type detachable, 5 Lead
E202-5029/XX	Dual pole detachable, 5 Lead	
 Connector No: 037 MENNEN 900, 1100, 1110, Horizon 2000, 2200, 5200, 9000XL, MR1200, 1270, 1280, 1300, 1380	E200-3037/XX	One piece (monolithic), 3 Lead
	E201-3037/XX	DIN type detachable, 3 Lead
	E202-3037/XX	Dual pole detachable, 3 Lead
	E200-5037/XX	One piece (monolithic), 5 Lead
	E201-5037/XX	DIN type detachable, 5 Lead
E202-5037/XX	Dual pole detachable, 5 Lead	
 Connector No: 075 MENNEN Envoy, Vitalogic 4000/5000, Vitalogic 5000/5500 – configured monitor	E200-3075/XX	One piece (monolithic), 3 Lead
	E201-3075/XX	DIN type detachable, 3 Lead
	E202-3075/XX	Dual pole detachable, 3 Lead
	E200-5075/XX	One piece (monolithic), 5 Lead
	E201-5075/XX	DIN type detachable, 5 Lead
E202-5075/XX	Dual pole detachable, 5 Lead	



Table-2 ECG Patient Cables Compatibility Table (Continued)

Connector	Part Number	Yoke Type
 Connector No: 015  METRAX PRIMEDIC Defi N/B, Defimonitor ECO1, DM1/3/10/30	E200-4015/XX	One piece (monolithic), 4 Lead
	E201-4015/XX	DIN type detachable, 4 Lead
	E202-4015/XX	Dual pole detachable, 4 Lead
 Connector No: 065  METRAX PRIMEDIC	E200-4065/XX	One piece (monolithic), 4 Lead
	E201-4065/XX	DIN type detachable, 4 Lead
	E202-4065/XX	Dual pole detachable, 4 Lead
 Connector No: 033  METRAX PRIMEDIC Hardsave 6/6S, Defimonitor XD.	E200-3033/XX	One piece (monolithic), 3 Lead
 Connector No: 169  METRAX PRIMEDIC Defimonitor XD, XD1, XD3, XD300.	E200-4016/XX	One piece (monolithic), 4 Lead
 Connector No: 001  METSIS Life Point (Monophasic)	E200-5001/XX	One piece (monolithic), 5 Lead
	E201-5001/XX	DIN type detachable, 5 Lead
	E202-5001/XX	Dual pole detachable, 5 Lead
 Connector No: 133  METSIS Life Point	E200-5071/XX	One piece (monolithic), 5 Lead
	E201-5071/XX	DIN type detachable, 5 Lead
	E202-5071/XX	Dual pole detachable, 5 Lead
 Connector No: 098  METSIS Life Point (Biphasic)	E200-4027/XX	One piece (monolithic), 4 Lead
	E201-4027/XX	DIN type detachable, 4 Lead
	E202-4027/XX	Dual pole detachable, 4 Lead
 Connector No: 001  MINDRAY PM9000, PM7000, PM8000	E200-3001/XX	One piece (monolithic), 3 Lead
	E201-3001/XX	DIN type detachable, 3 Lead
	E202-3001/XX	Dual pole detachable, 3 Lead
	E200-5001/XX	One piece (monolithic), 5 Lead
	E201-5001/XX	DIN type detachable, 5 Lead
	E202-5001/XX	Dual pole detachable, 5 Lead



Table-2 ECG Patient Cables Compatibility Table (Continued)

Connector	Part Number	Yoke Type
 Connector No: 084 MINDRAY	E200-3084/XX	One piece (monolithic), 3 Lead
	E201-3084/XX	DIN type detachable, 3 Lead
	E202-3084/XX	Dual pole detachable, 3 Lead
	E200-5084/XX	One piece (monolithic), 5 Lead
	E201-5084/XX	DIN type detachable, 5 Lead
	E202-5084/XX	Dual pole detachable, 5 Lead
PM5000, PM6000, PM8000.		
 Connector No: 245 NEC	E200-3045/XX	One piece (monolithic), 3 Lead
	E201-3045/XX	DIN type detachable, 3 Lead
	E202-3045/XX	Dual pole detachable, 3 Lead
	E200-5045/XX	One piece (monolithic), 5 Lead
	E201-5045/XX	DIN type detachable, 5 Lead
	E202-5045/XX	Dual pole detachable, 5 Lead
 Connector No: 009 NIHON KOHDEN	E200-3009/XX	One piece (monolithic), 3 Lead
	E201-3009/XX	DIN type detachable, 3 Lead
	E202-3009/XX	Dual pole detachable, 3 Lead
	E200-5009/XX	One piece (monolithic), 5 Lead
	E201-5009/XX	DIN type detachable, 5 Lead
	E202-5009/XX	Dual pole detachable, 5 Lead
TEC-5500/5521/5531/7600/7700.		
 Connector No: 099 NIHON KOHDEN	E200-3099/XX	One piece (monolithic), 3 Lead
	E201-3099/XX	DIN type detachable, 3 Lead
	E202-3099/XX	Dual pole detachable, 3 Lead
 Connector No: 090 NIHON KOHDEN	E200-3090/XX	One piece (monolithic), 3 Lead
	E201-3090/XX	DIN type detachable, 3 Lead
	E202-3090/XX	Dual pole detachable, 3 Lead
	E207-3090/XX	Eight pole Nihon Kohden, 3 Lead
	E207-6090/XX	Eight pole Nihon Kohden, 6 Lead
	BSM 2300/2301/4100/5100/9510/9800, Lifescope S, Lifescope N.	
 Connector No: 204 NIHON KOHDEN	E200-3086/XX	One piece (monolithic), 3 Lead
	E201-3086/XX	DIN type detachable, 3 Lead
	E202-3086/XX	Dual pole detachable, 3 Lead
	E207-3086/XX	Eight pole Nihon Kohden, 3 Lead
	OPV 1500	
 Connector No: 008 NIHON KOHDEN	E200-3008/XX	One piece (monolithic), 3 Lead
	E201-3008/XX	DIN type detachable, 3 Lead
	E202-3008/XX	Dual pole detachable, 3 Lead
	E200-5008/XX	One piece (monolithic), 5 Lead
	E201-5008/XX	DIN type detachable, 5 Lead
	E202-5008/XX	Dual pole detachable, 5 Lead
OEC-6102A(S/N<20156), OEC-6105A(S/N<20226), OEC-6120A/J/K		
 Connector No: 001 NIHON KOHDEN	E200-3001/XX	One piece (monolithic), 3 Lead
	E201-3001/XX	DIN type detachable, 3 Lead
	E202-3001/XX	Dual pole detachable, 3 Lead
	E200-5001/XX	One piece (monolithic), 5 Lead
	E201-5001/XX	DIN type detachable, 5 Lead
	E202-5001/XX	Dual pole detachable, 5 Lead
BSM-8301A, 8302, 8500A, 8502, OEC6102A, 6105A, 7102A, 8108A, TEC 5200A, 7100A, 7200A, 7300A, 7541A		



Table-2 ECG Patient Cables Compatibility Table (Continued)









Connector	Part Number	Yoke Type	
 Connector No: 001 PACE-TECH	 Minipack 300/3000/3100, Vitalmax 800/2100/2200/4000/4100.	E200-3001/XX	One piece (monolithic), 3 Lead
		E201-3001/XX	DIN type detachable, 3 Lead
		E202-3001/XX	Dual pole detachable, 3 Lead
		E200-5001/XX	One piece (monolithic), 5 Lead
		E201-5001/XX	DIN type detachable, 5 Lead
		E202-5001/XX	Dual pole detachable, 5 Lead
 Connector No: 004 PALCO	 500, 500VP.	E200-5044/XX	One piece (monolithic), 5 Lead
		E201-5044/XX	DIN type detachable, 5 Lead
		E202-5044/XX	Dual pole detachable, 5 Lead
 Connector No: 140 PETAŞ	 PM100, KMA275/375/675/800.	E200-3003/XX	One piece (monolithic), 3 Lead
		E201-3003/XX	DIN type detachable, 3 Lead
		E202-3003/XX	Dual pole detachable, 3 Lead
		E200-4003/XX	One piece (monolithic), 4 Lead
		E200-5003/XX	One piece (monolithic), 5 Lead
		E201-5003/XX	DIN type detachable, 5 Lead
 Connector No: 004 PETAŞ	 KMA 165/265/365.	E200-4004/XX	One piece (monolithic), 3 Lead
		E201-4004/XX	DIN type detachable, 3 Lead
		E202-4004/XX	Dual pole detachable, 3 Lead
		E200-5004/XX	One piece (monolithic), 5 Lead
		E201-5004/XX	DIN type detachable, 5 Lead
		E202-5004/XX	Dual pole detachable, 5 Lead
 Connector No: 020 PHILIPS / HP	 Please see Note 1	E200-3020/XX	One piece (monolithic), 3 Lead
		E201-3020/XX	DIN type detachable, 3 Lead
		E202-3020/XX	Dual pole detachable, 3 Lead
		E205-3020/XX	Dual pole detachable, 3 Lead
		E200-5020/XX	One piece (monolithic), 5 Lead
		E201-5020/XX	DIN type detachable, 5 Lead
 Connector No: 024 PHILIPS / HP	 Please see Note 2	E200-3024/XX	One piece (monolithic), 3 Lead
		E201-3024/XX	DIN type detachable, 3 Lead
		E202-3024/XX	Dual pole detachable, 3 Lead
		E205-3024/XX	Dual pole detachable, 3 Lead
		E200-5024/XX	One piece (monolithic), 5 Lead
		E201-5024/XX	DIN type detachable, 5 Lead
 Connector No: 022 PROGETTI	 EPG View Defibrillator	E200-3035/XX	One piece (monolithic), 3 Lead
		E200-5035/XX	One piece (monolithic), 5 Lead

Note 1 ;
 43100A/43110A/43110M/43120A/MC, 43200A, M3500B, 7700AC/77020A/ 77010AC, 77025A, 77030A, 77035A, 77040A, 78208A, 78210A, 7330A, 78331A, 78332A, 78333A, 78341A, 78342A, 78345A, 78341A, 78351T, 78352A/C, 78353A/B, 78354A, 78532A/B, 78354A/B/C, 78619A, 78620A, 78660A/B, 78661A/B, 78670A, 78671A, M1722A/B, M2400A, M2406A, Donor 2000/2500, M2425A, M2475B

Note 2 ;
 78352C, 78354C, 78834C, 862474, C3, M1001A, M1001A, M1002A, M1165A, M1166A, M11767A, M1175A, M1177A, M1205A, M1264A, M1275A, M1280A, M1723A/3, M2475B, M2601A, M3000A, M30001A, M3046A, M4M3, M4735A



Table-2 ECG Patient Cables Compatibility Table (Continued)

Connector	Part Number	Yoke Type
 Connector No: 022 RGB Omicron FT Surveyor, Omicron Vision	E200-3092/XX	One piece (monolithic), 3 Lead
	E201-3092/XX	DIN type detachable, 3 Lead
	E202-3092/XX	Dual pole detachable, 3 Lead
	E200-5092/XX	One piece (monolithic), 5 Lead
	E201-5092/XX	DIN type detachable, 5 Lead
	E202-5092/XX	Dual pole detachable, 5 Lead
 Connector No: 085 SAADAT Novin S1600, Novin S1800.	E200-3085/XX	One piece (monolithic), 3 Lead
	E201-3085/XX	DIN type detachable, 3 Lead
	E202-3085/XX	Dual pole detachable, 3 Lead
	E200-5085/XX	One piece (monolithic), 5 Lead
	E201-5085/XX	DIN type detachable, 5 Lead
	E202-5085/XX	Dual pole detachable, 5 Lead
 Connector No: 022 SCHILLER Argos LCM.	E200-3081/XX	One piece (monolithic), 3 Lead
	E201-3081/XX	DIN type detachable, 3 Lead
	E202-3081/XX	Dual pole detachable, 3 Lead
	E200-4081/XX	One piece (monolithic), 4 Lead
	E201-4081/XX	DIN type detachable, 4 Lead
	E202-4081/XX	Dual pole detachable, 4 Lead
	E200-5081/XX	One piece (monolithic), 5 Lead
	E201-5081/XX	DIN type detachable, 5 Lead
	E202-5081/XX	Dual pole detachable, 5 Lead
 Connector No: 022 SCHILLER Defiguard 2002/3002/4000/5000.	E200-3082/XX	One piece (monolithic), 3 Lead
	E200-4082/XX	One piece (monolithic), 4 Lead
 Connector No: 030 SIEMENS - DRÄGER Sinecast, 400, 500, 600, 700, 900, 1200, 341, 610, 620, 630, 720, 722, 730, 732, 960, 961, 1260, 1281, 1801.	E200-3030/XX	One piece (monolithic), 3 Lead
	E201-3030/XX	DIN type detachable, 3 Lead
	E202-3030/XX	Dual pole detachable, 3 Lead
	E200-5030/XX	One piece (monolithic), 5 Lead
	E201-5030/XX	DIN type detachable, 5 Lead
	E202-5030/XX	Dual pole detachable, 5 Lead
 Connector No: 083 SIEMENS - DRÄGER Drager VM1/3, RM1, PM8010, PM8014, Babyguard 8000, Cato, Babylog 8000/8000 plus, PB8800, Cicero.	E200-3083/XX	One piece (monolithic), 3 Lead
	E201-4083/XX	DIN type detachable, 3 Lead
	E202-5083/XX	Dual pole detachable, 3 Lead
 Connector No: 038 SPACELABS Ultraview 90496, 90369, 90367.	E200-3038/XX	One piece (monolithic), 3 Lead
	E201-3038/XX	DIN type detachable, 3 Lead
	E202-3038/XX	Dual pole detachable, 3 Lead
	E206-3038/XX	One piece (monolithic), 3 Lead
	E200-5038/XX	One piece (monolithic), 5 Lead
	E201-5038/XX	DIN type detachable, 5 Lead
	E202-5038/XX	Dual pole detachable, 5 Lead
	E206-5038/XX	One piece (monolithic), 5 Lead
 Connector No: 001 SPACELABS Please See Note 3	E200-3001/XX	One piece (monolithic), 3 Lead
	E201-3001/XX	DIN type detachable, 3 Lead
	E202-3001/XX	Dual pole detachable, 3 Lead
	E200-5001/XX	One piece (monolithic), 5 Lead
	E201-5001/XX	DIN type detachable, 5 Lead
	E202-5001/XX	Dual pole detachable, 5 Lead

Note 3 :
 400 Series (5 Id), 900 Series (w/respiration), 90300 / 90303B / 90305 / 90308 / 90408 / 90418 / 90425 / 90428 / 90432 / 90470, 90700 Series, TEK-208/408/413/414 (w/lead set), 500/600 Series, 4045 Defib, PC Bedside, PCMS System.



Table-2 ECG Patient Cables Compatibility Table (Continued)

Connector	Part Number	Yoke Type
 <p>Connector No: 153</p> <p>SPENGLER</p> <p><i>All Models.</i></p>	E200-3070/XX	One piece (monolithic), 3 Lead
	E201-3070/XX	DIN type detachable, 3 Lead
	E202-3070/XX	Dual pole detachable, 3 Lead
	E200-5070/XX	One piece (monolithic), 5 Lead
	E201-5070/XX	DIN type detachable, 5 Lead
	E202-5070/XX	Dual pole detachable, 5 Lead
 <p>Connector No: 010</p> <p>S&W (ARTEMA)</p> <p><i>Diascope 1/2, Diascope NT, Athena, Cardio-Aid, Traveller, Neoscope.</i></p>	E200-3010/XX	One piece (monolithic), 3 Lead
	E201-3010/XX	DIN type detachable, 3 Lead
	E202-3010/XX	Dual pole detachable, 3 Lead
	E200-5010/XX	One piece (monolithic), 5 Lead
	E201-5010/XX	DIN type detachable, 5 Lead
	E202-5010/XX	Dual pole detachable, 5 Lead
 <p>Connector No: 096</p> <p>UTAS</p>	E200-5096/XX	One piece (monolithic), 5 Lead
	E201-5096/XX	DIN type detachable, 5 Lead
	E202-5096/XX	Dual pole detachable, 5 Lead
 <p>Connector No: 001</p> <p>WELLCH ALLYN</p> <p><i>Propaq 100/200 Series, Encore, CS; Ultra Smart.</i></p>	E200-3002/XX	One piece (monolithic), 3 Lead
	E201-3002/XX	DIN type detachable, 3 Lead
	E202-3002/XX	Dual pole detachable, 3 Lead
	E200-5002/XX	One piece (monolithic), 5 Lead
	E201-5002/XX	DIN type detachable, 5 Lead
	E202-5002/XX	Dual pole detachable, 5 Lead
 <p>Connector No: 001</p> <p>ZOLL</p> <p><i>D900, PD1200/1400/1600</i></p>	E200-3002/XX	One piece (monolithic), 3 Lead
	E200-5002/XX	One piece (monolithic), 5 Lead

Multi Parameter Cables

FMT-MTM/01 (OEM P/N # 3368391) Multi Parameter Cable is designed to monitor 3 or 5 lead ECG, impedance respiration, temperature and SpO2.

FMT-MTM/01 Multi Parameter Cables are compatible with SC6000, SC7000, SC8000, SC9000, Drager Infinity, Delta and Delta XL monitors (with Nellcor SpO2 option).



Holter ECG Leadwires

HL01-7

Holter leadwire, 7 leadwire set, 1.5 mm safety DIN, snap patient end, 90 cm, AHA

HL01-10

Holter leadwire, 10 leadwire set, 1.5 mm safety DIN, snap patient end, 90 cm, AHA



HL02


Mortara holter leadwire, 10 leadwire set, snap patient end, IEC



Prewired Neonatal/Pediatric ECG Electrodes

FMT NE Neonatal ECG Electrodes are Ø22 mm in diameter and consist of a conductive adhesive hydro gel, a Silver-Silver Chloride (Ag/AgCl) plated patient contact, a non-woven substrate and a 60cm flexible tinned copper or carbon lead-wire terminating in a standard Safety 1.5 mm DIN style connector.

NE Neonatal ECG Electrodes supplied 3 pieces set in a pouch. Noise free traces are obtained with the high quality very low impedance, hypoallergenic gel. Radiolucent version of the prewired neonatal/pediatric ECG electrodes is also available upon request.




NEI (IEC Green Yellow Red), **NEA** (AHA Red, Green, Black), **NEW** (All White), **NET** (Carbon conductor, Transparent)




FMT

HASTABAŞI MONİTÖRÜ EKG KABLOLARI ÜRÜN DÖKÜMANI

► Renk Kodları ve Dizilimi



FMT® Metko Ltd.'nin tescilli ticari markasıdır. Burada bahsedilen diğer ürünler, logolar ve şirket isimleri başka firmaların ticari markaları olabilir. Bu durum yukarıda geçen markaların FMT® ile herhangi bir yakınlığı veya onayı olarak algılanmamalıdır. Bu dokümanda verilen bilgiler dikkatlice kontrol edilmiş ve doğruluğuna inanılmıştır. Sürekli ürün gelişimi adına, Metko bu dokümanda ve tarif edilen üründe uyarı ve bilgilendirme yapmadan değişiklikler ve gelişmeler yapma hakkını saklı tutar.



HEADQUARTER & FACTORY

- ▶ İvedik O.S.B. Ağaç İşleri Sanayi Sit.
1354. Cad. 1358.Sok. No:9
06378 Yenimahalle

ANKARA / TÜRKİYE

Tel (pbx) : +90 312 387 12 46
Fax : +90 312 387 12 51

BRANCH

- ▶ Atatürk Mah. Sedef Cad. 38 Ada
Ata Plaza 3-5 Blokları Kat:2 No:57
34758 Ataşehir

İSTANBUL / TÜRKİYE

Tel (pbx) : +90 216 455 91 21
Fax : +90 216 455 48 24

INTERNET

- ▶ <http://www.metkoltd.com>
metko@metkoltd.com